

Cervical Cancer Screening Sample Collection

FJORD HPV-SAMPLER™

Designed for sampling, transportation, and storage of HPV DNA test sample

1. FJORD HPV-SAMPLER™ - Cervical Sample

Specimen type: Cervical
(Endocervical or Ectocervical)



2. FJORD HPV-SAMPLER™ - Urine Sample

Specimen type: Urine
A preferred sample type of collection for women

(Urine HPV viral load is processed with SEDIPREP® Pre-Concentration prior to molecular analysis.)



A Proven Choice for Cervical & Urine Specimen Storage, Transportation and Processing

Cervical Cancer Screening Sample Collection

Product Details

FJORD HPV-SAMPLER™ is a liquid-based cytology collection solution that provides long term stability for specimen containing cellular materials with human papillomavirus, to be used in molecular HPV DNA testing and/or for the preparation of cytology slides, enabling the support for the full spectrum of cervical cancer screening workflow.

Performance

Specimen integrity is retained for 21 days at room temperature with minimal degradation.

Compatibility

Validated for use with molecular assays, like HC2 HPV DNA, Cobas 4800 HPV, etc.

Screening Applications

	FJORD HPV-SAMPLER™	
	<i>Cervical Sample</i>	<i>Urine Sample</i>
ASCUS triage	✓	
Co-testing (Pap Cytology with HPV testing)	✓	
Primary HPV Screening	✓	✓
Pap Cytology - CY-PREP® Pap Test (or compatible LBC Pap Test)	✓	