


THE CY-PREP® Dual Filtration Liquid Based Cytology System

- ✓ Improving diagnostics for clinicians
- ✓ Improving profitability for laboratories

THE CY-PREP® ADVANTAGE

Unwanted debris in liquid cytology specimen are often the source of interfering elements in the preparation of consistent quality slides.

The CY-PREP® Dual Filtration Process effectively removed these unwanted debris at the 1st filtration membrane process, leaving only the wanted cell materials that are captured and retained on the surface of the 2nd filtration membrane for transfer onto a glass slide in a thin layer format.



The CY-PREP® system automates every step of the slide preparation process. Slides are consistently prepared with excellent cellularity, morphology and distribution.